# T° PATCH

## PACKAGE LEAFLET

ITALIANO

ENGLISH

DEUTSCH

FRANÇAIS

ESPAÑOL

PORTUGUÊS

### **Descrizione:**

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## PACKAGE LEAFLET

### Description:

**T° Patch** is a Class I EC medical device, in accordance with Annex VII of MDD 93/42/EEC. Pioneering the new concept of the wearable thermometer, the T°Patch is an adhesive thermochromic temperature indicator that is as easy to apply as a plaster, conformable to the skin and water resistant. Its hi-tech thermochromic paint allows an immediate and constant check of body temperature. Its three colour codes make the device easy to use regardless of the age and health condition of the consumer. Its water-resistant materials allow prolonged use of the device for up to 4 days for adults and up to 2 days for children. For an accurate indication of body temperature, we recommend placement on the neck, chest, underarm or forehead.

### Technical characteristics:

- Class I EC medical device
- Non-sterile and non-reusable product
- Accuracy of  $\pm 0.2$  °C.
- Waterproof
- Recommended period of use: up to 4 days for adults and up to 2 days for children.

### Color Code:

Figure 1 Color code

	Color	Temperature
LOW	Brown	$T < 34\text{ }^{\circ}\text{C}$
NORMAL	Green	$34\text{ }^{\circ}\text{C} < T < 37\text{ }^{\circ}\text{C}$
HIGH	Yellow	$T > 37.5\text{ }^{\circ}\text{C}$

The colour of the T° Patch varies according to body temperature with an accuracy of  $\pm 0.2\text{ }^{\circ}\text{C}$ . The T° Patch will be **brown** with a body temperature up to  $34\text{ }^{\circ}\text{C}$ , **green** with a temperature between  $34\text{ }^{\circ}\text{C}$  and  $37\text{ }^{\circ}\text{C}$  and **bright yellow** with a temperature above  $37.5\text{ }^{\circ}\text{C}$ .

#### Use:

1. Open the package.
2. Remove one of the eight T°Patches and apply it to one of the recommended areas (see section "Application areas").
3. Check the colour of the T°Patch during the day.
4. Remove the T°Patch after the fourth day for adults, after two days for children.

## Application Areas:

Scientific studies show that human body temperature is not constant in all parts of the body, but varies according to the flow of blood to that particular body area (Fig. 2).



Figure 2 Termochromic map of the human body

Based on these studies (see bibliography), the areas recommended for the application of the T° Patch are as follows:

1. Neck.
2. Behind the ear.
3. Breast.
4. Forehead.



Figure 3 Recommended application areas

## Warnings:

- Apply after thoroughly cleaning the skin.
- Do not apply to areas covered with hair.
- Do not ingest.
- Keep out of reach of children under 36 months.
- Do not use after the expiry date.
- For temperature quantification use a standard thermometer.
- Non-sterile product.
- Not reusable product.

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[2] Senkevich, Y.I., Yafarov, A.Z. Simulation of Human Body Temperature Variation under Cold Conditions. *Biomed Eng* **48**, 93–96 (2014). <https://doi.org/10.1007/s10527-014-9428-z>

[3] Boratyński, J.S., Iwińska, K. & Bogdanowicz, W. Body temperature variation in free-living and food-deprived yellow-necked mice sustains an adaptive framework for endothermic thermoregulation. *Mamm Res* **63**, 493–500 (2018). <https://doi.org/10.1007/s13364-018-0392-y>

[4] Levin, Y.A., Prikhod'ko, I.M. Temperature of a thin body with periodic variation of the heat transfer coefficient and the temperature of the surrounding medium. *Journal of Engineering Physics* **10**, 144–145 (1966). <https://doi.org/10.1007/BF00834556>